


# Hearing Aids

---

Carol De Filippo

Viet Nam Teacher Education Institute

June 2010





# Topics: Hearing Aids

- Benefits
- Components
- Technologies
- Styles
- Fitting
- Care
- Classroom teacher tasks

# Hearing Aids

... for *all* degrees of hearing loss



- **Mild Hearing Loss**

- Very helpful for hearing high-frequency speech sounds (examples: p, t, k, f, s, sh)
- Especially beneficial when hearing loss reaches 30-40 dB HL

- **Moderate to Severe Hearing Loss**

- High motivation
- High benefit



## ● **Profound Hearing Loss**

- Great variation in benefit
  - May shift functional status to “hard-of-hearing”
  - May not be a choice for those who identify as culturally Deaf
- Factors that influence benefit
  - Age when hearing loss was identified
  - Type of hearing loss and progression of loss
  - Age when hearing aid was first used (duration of deafness)
  - Exposure to spoken language
  - Consistency of listening experience
  - Family support
  - School support



# Benefits of Hearing Aids

- Expanded **communication** options
  - Greater hearing aid use correlates with increased speech intelligibility
  - Greater hearing aid use correlates with more effective use of speechreading
- **Education** gains
  - Continued hearing aid use correlates with higher academic achievement
  - Duration of hearing aid use correlates with language skill level




<http://www.time4speech.com.au/Page2-Speech-Sounds/Speech-Sounds.html>

# Reasons for nonuse in the U.S.



- Lack of encouragement
- Lack of models (successful adults in the environment with hearing aids)
- Lack of support
- Lack of auditory environment (no *need to listen*)
- Older children's concern that a hearing is very noticeable; peer pressure
- Poor hearing aid fit, which leads to poor benefit

# Basic hearing aid components





- **Microphone**

- Collects sound waves from the air
- Convert waves to an electronic signal

- **Signal processor**


- Increases amplitude of the signal
- Modifies signal to ensure audibility, comfort, and protection from loud sound

- **Receiver**

- Converts electronic signal back to acoustic sound waves

- **Battery**

- Provides power





# Hearing aid technology




- **Digital signal processing**
  - Precise control of all fitting parameters
  - Takes up less space, so aid can be smaller
  - Analyzes incoming signals and then selects the appropriate processing strategy
    - Examples: automatic noise reduction, automatic speech compression
- **Multiple programs**
  - Hearing aid settings are stored for several different listening situations
  - Examples: one-to-one quiet conversation, large-group meeting, telephone
- **Bluetooth connectivity**
  - Wireless communication between a hearing aid and another device (examples: cell phone, music player, hearing aid on the other ear)
- **Remote control** to select programs and settings




# Hearing aid technology

- **Microphones**

- Omni-directional: Takes in sound from all directions
- Directional: Takes in sound from all directions
- Demonstration



Car noise interferes




Hearing aid focuses on the person's speech

# Hearing aid technology




- **Compression amplification**

- The hearing aid automatically fits the speech signal into the listener's range of hearing
  - Above the threshold of detection
  - Below the threshold of discomfort
  - Without distorting the quality of the sound



- Amount of amplification adapts to the level of the input

- Weak inputs are amplified more
- High-level inputs are amplified less




# Hearing aid technology



- **Feedback cancellation**

- A squealing sound made by the hearing aid
- Caused when the amplified sound is re-amplified
- Can be controlled by digital sound processing



- **Telephone connection**

- Telecoil
  - A special circuit that picks up and amplifies the electromagnetic signals that come out of a telephone handset (or any loop of wire)
  - Blocks out surrounding noise







# Hearing aid technology

- **Direct audio input**

- To connect to devices such as TV, telephone, computer, CD player
- Hearing aid has an external connector or contacts
  - To insert a mini-plug to another device
  - To attach a “boot” and mini-plug that connects to another device
  - To attach a “boot” with a tiny FM radio receiver that picks up the signal from another device
- Improves the signal because it bypasses the hearing aid microphone and eliminates background noise



# Hearing aid technology



- **Earmold** materials and features

- A good earmold reduces many hearing aid problems

- **Materials**

- Silicone
  - Softest
  - For sports; to control acoustic feedback (squeal)
- Lucite
  - Hard; easily modified
  - More durable
  - NOT used with young children to avoid injury if struck on the side of the head


- **Earmold tubing**

- Shape and thickness can change the sound output
- Will harden
  - Needs to be changed every 6-12 months

- **Earmold venting**



- Reduces occlusion effect (annoyance from talking, chewing)





# Styles of Hearing Aids


- In-the-ear (ITE)
- Behind-the-ear (BTE)
- Body aid
- Bone-Anchored Hearing Aid (BAHA)
- Tactile aids



# In-the-Ear Hearing Aids


## Advantages

- Small size
- Comfort (depends on the individual)
- Not very noticeable
- Sound quality
- Very little wind noise
- CIC & ITC can reduce occlusion effect




Up to 50 dB

Completely-in-the-canal (CIC)



Up to 50 dB

In-the-canal (ITC)



Up to 70 dB

In-the-ear (ITE)





# In-the-Ear Hearing Aids

## *Disadvantages*

- Cannot fit severe or profound losses
- More expensive (custom fit housing)
- More difficult to manipulate
- May have fewer features (no telephone connection, directional microphone, venting)
- May need remote control device
- Prone to wax blockage
- May be uncomfortable
- Easier to lose
- Need loaner when out for repair



CIC



ITC



ITE



# Behind-the-Ear Hearing Aids



## *Advantages*

- High power for profound hearing loss
- Least expensive
- Most flexible to fit a wide range of needs
- User can wear a loaner aid with his own earmold when repairs are needed
- Multiple features
  - Direct audio input or FM input
  - Telephone connection (telecoil)
  - Program switch
  - Directional microphone
  - Durable
  - Easy to manipulate
  - Longest battery life
  - Full earmold venting possible





## *Disadvantages*

- Less comfortable on the ear (depends on the individual)
- Wind noise (depends on location of the microphone)
- More noticeable





# Advances in BTE design

## ● Thin-tube, open-ear fit

- Avoids a fitted earmold and leaves the ear canal more open
- The more open the mold, the more comfortable
- Open system advantages:
  - Ventilation of the ear canal
  - Avoids overheating
  - Avoids accumulation of moisture
  - Creates more natural sounding voice
  - Potential for optimal amplification in the high frequencies
- Works for mild & moderate losses



CONVENTIONAL

OPEN-EAR FIT



# Advances in BTE design

- **Receiver-in-the ear (RITE)**
  - Avoids acoustic resonance of earmold tubing
  - Places sound further down the ear canal
  - Can improve clarity of the signal
    - Can provide a smooth frequency response



[www.oticonusa.com/.../HearingAidSelection.html](http://www.oticonusa.com/.../HearingAidSelection.html)



Closed dome for greater hearing loss



1. Amplifier
2. Thin sound wire
3. Dome
4. Receiver
5. Ear grip

# Body Aids



## *Advantages*

- Highest gain (example: up to 130 dB)
- Controls are easy to manipulate (for small children; elderly people)
- Has many features with flexible fitting options

## *Disadvantages*

- Cords can break
- Cord and “box” are very noticeable
- No directional microphone
- May not be programmable
- May not have a digital processor






# Tactile Aids

## *Design*

- Vibrators are used to deliver the acoustic signal to the skin (fingers, forearm, torso, or back)
- Example: Two vibrators might represent low frequencies (vowels) and high frequencies (consonants)

## *Advantages*

- Provides access to the rhythmic pattern of speech
- Beneficial as a lipreading aid
- Helpful for speech training



## *Disadvantages*

- Useful only in quiet settings
- Cannot obtain detailed frequency information

# Bone-Anchored Hearing Aid

## BAHA



- For persons with conductive or mixed hearing loss (outer or middle ear problems)
  - Chronic otitis media (middle ear infection)
  - Microtia or other malformations of the outer ear or middle ear (congenital or due to trauma or cancer)
- A surgically implanted system




Abutment



Processor

<http://www.umm.edu/otolaryngology/baha.htm>



<http://products.cochlearamericas.com/baha/introduction-to-baha>



# Binaural hearing aids



## ● The binaural advantage

- Eliminates head shadow
  - The listener can receive sound coming from either side at a comfortable level
- Loudness summation
  - A listener with high thresholds can get a "boost" by using two aids
  - Less amplification is needed because a signal received by two ears is heard as "louder" than the same signal received by one ear
- Binaural squelch
  - It is much easier to pick a signal out of the noise when the brain receives information from both ears
- Localization
  - The brain can calculate the direction of the signal with information from two ears (localization is impossible with input from one ear only)
- Convenience
  - If one hearing aid needs repair, the listener does not have to go without amplification



- **Disadvantages**

- All costs are doubled
  - Purchase costs
  - Maintenance and repair costs
- The combined signals may not blend
  - There may be a significant difference between ears in threshold or sound quality
  - The brain may require a lot of training to "pay attention" to the new input if one ear has not been stimulated in a long time
  - Neural capacity to use the input may have diminished

# Hearing aid fitting



## 1. Fitting Stage

- Select a hearing aid with the desired gain and other features
- Use a computerized program to adjust the hearing aid to a “prescription” (rules for setting hearing aid levels)
  - Provide enough amplification so that soft sounds can be heard
  - Provide limits so that loud sounds are not uncomfortable

## 2. Verification Stage

- Measure electroacoustic properties of the hearing aid
- Perform “real-ear measurements”
  - Place small microphone in ear canal
  - Measure sound reception in open ear canal
  - Measure sound reception in ear canal with HA in place
  - Calculate difference with and without HA
  - Compare gain to amplification prescription
  - Adjust the HA as needed to match prescription





### 3. **Validation** stage

- Measure benefit
  - Speech recognition tests
  - Parent report
  - Self report (questionnaires, inventories)



# Hearing aid care

- **Cautions**

- Do not place a hearing aid in direct heat or sunlight (example, hot car)
- Avoid dropping the aid
- Avoid getting the aid wet (no swimming or showering)
- Keep aid away from pets and small children

- **Daily care** of the hearing aid

- Wipe off any dirt or moisture with a soft, dry cloth
- Remove earwax with special wire pick
- Remove battery or open battery compartment overnight
- Store overnight in a “dry-sac” (with a desiccant to absorb moisture)
- Check battery voltage and replace battery when low
  - Remove dead batteries promptly
  - Turn aid off when not in use

# Hearing aid care



- **BTE earmold and tubing**

- Life expectancy of earmold/tubing is 1 year for an adult, 3 months for a child
- Material shrinks over time, which will cause feedback
- Clean every 1-7 days
  - Use *warm* water & soap
  - Remove all moisture by blowing air through the tubing



# Classroom hearing aid checks



## 1. Visual Inspection

- Check earmold
  - Moisture
  - Wax
- Inspect the case
  - Cracks
  - Broken or missing parts
- Inspect hook and earmold tubing
  - Hard
  - Cracked, split, punctured
  - Twisted tubing
- Check battery and battery compartment
  - Match + to +
  - Measure voltage
  - Clean any corrosion with sandpaper or eraser
- Inspect cords for visible breaks

## 2. Listening Check #1

- Turn hearing aid volume full on
- Plug end of earmold
- Check for sound leakage

## 3. Listening Check #2

- Attach a stethoset
- Listen again
- Carefully increase and decrease volume
  - Listen for breaks, buzz, crackling noises



# Classroom hearing aid checks



## 4. Check Student Function

- Use the Ling 7-sound test
  - **oo, ah, ee, sh, m, s, or** (example for English)
  - These sounds represent the variety of frequencies present in speech
- Set a baseline for each student
  - Use a functioning hearing aid with a good battery
  - Record which sounds the student can perceive from 3 feet away
  - Not all students will hear all 7 sounds
- Daily check
  - Student sits or stands about 3 feet away, wearing his/her hearing aid
  - Cover your mouth with a card
  - Individually present each of the 7 sounds
  - Have the student raise a hand, drop a block into a container, or say "yes," if the sound is heard





# Review: Hearing Aids

- Benefits
- Components
- Technologies
- Styles
- Fitting
- Care
- Classroom teacher tasks