

## Topspin Protocol for Bruker 500 NMR

### I. Sample Introduction

- a. Login to the NMR workstation using the proper username and password
- b. On the desktop open “Topspin 3.0”.
- c. Type “**ej**” in the bottom window of Topspin.
  - i. Note: Wait for the air to turn on before inserting sample.
- d. Insert the sample into the magnet
- e. Type “**ij**” to lower sample into the magnet
- f. Type “**ro**” to spin sample

### II. Experimental Setup

- a. At the bottom of topspin type “**rsh latest**”
- b. Click **File** then **New** on the top menu.
  - i. **Name** – Define a file name for your sample
  - ii. **Expno** – This is the experimental number and should begin with **1**
  - iii. **Procno** – this is the process number, normally set to **1**
  - iv. **Dir** – This should be set to C:\Bruker\Topspin 3.0\Data
  - v. **User** – This should be the same username you are logged in under
  - vi. **Solvent** – Select the deuterated solvent used.
  - vii. **Experiment** – Select **Proton** from the dropdown menu.
  - viii. **Title** – Enter a title for the spectrum. This will be displayed in the printout.
- c. Click on **OK** in the window to close

### III. Locking, Shimming, and Wobbling

- a. Type “**lock**” and choose the appropriate lock solvent
  - i. When lock is finished, as noted in the bottom window, proceed to the next step.
- b. Type “**getprosol**” in the bottom window of Topspin
- c. Type “**atma**”, in the bottom window to wobb the sample.
- d. When wobbling is finished, as indicated at the bottom of topspin, then type “**topshim**” to shim the sample.
  - i. Note: Wait until shimming is completed before going on.


### IV. Data Acquisition







- a. Type “**rga**” and wait for finished to appear
- b. Type “**zg**”. You should see an FID signal after the dummy scans are completed

### V. Data Processing


- a. Once Data Acquisition is complete, type “**ef**” to perform exponential multiplication and a Fourier transform.
- b. Type “**apk**” to perform automatic phase correction.
- c. Type “**abs**” to perform automatic baseline correction.

### VI. Spectrum Calibration






- a. Click  on the top ribbon to begin calibration.

- i. Zoom in on the desired peak by adjusting:       or by clicking, holding, and dragging the cursor around the region desired.
- ii. Center the cursor on the peak to be calibrated, and click once.
- iii. Type the calibration frequency and click **OK**


#### VII. Automatic Peak Picking

- a. Click  then type “pps” to perform automatic peak picking.

#### VIII. Integrating Peaks

- a. Click  to begin integration.
  - i. Click  to select all integrals.
  - ii. Click  then **OK** to delete the automatically calculated integrals.
  - iii. Click  to begin defining new regions.
  - iv. Click, hold, drag, and release the cursor around the peak to be integrated. Repeat for remaining peaks.
  - v. Once finished with integration, click  to save changes.

#### IX. Plotting

- a. Type “**plot**” at the bottom of Topspin
- b. Double click on the center of the screen (green boxes should appear around the spectrum).
- c. Change Xmin / Xmax to the desired frequencies (in ppm).
  - i. Click **OK** when finished.
- d. Right mouse click on the center of the spectrum and select **1D/2D-Edit ...**
  - i. Adjust the vertical height using  and then click **Close** when finished.
- e. Click **File... Print** or press Ctrl+P to print.
- f. Close the Topspin Plot Editor and click **No** if a dialog asks if you wish to save changes.

#### X. Running a New Sample

- a. To run a new acquisition, start at step II, Experimental Setup

#### XI. To Eject a Sample

- a. To eject the sample type “**ej**” at the bottom of Topspin.
- b. Log off when finished